

Madison/REMD01

Electric rise dual tilt-in-space mechanism provides a range of postural positions, as the chair tilts and lifts the user into a standing position.



Standard Features

- Carer or user operated
- Separate waterfall back cushions can be individually repositioned and adjusted for comfort and support
- Seat depth adjustment - 1" or 2" adjuster pads
- Four heavy duty lockable castors
- Visco elastic pressure relieving foam seat
- Max load 125kg (20st)

Optional Features

- Additional pressure management seat cushions - Removable and interchangeable waterfall or lateral back cushions
- Alternative headrests and pillows
- Battery pack for Madison
- Arm covers, head cushion cover
- Honeycomb gel in leg & arm rests
- Lap strap with pad

DIMENSIONS



FINISH OPTIONS

- Available in a variety of healthcare fabric and vinyl options